

ATS™ Surgical Tables

ATS Trinity™

Carbon Fiber Top Imaging
& Vascular Table



Ideal for
vascular
centers

Intuitive Motion

Built to the foremost standard of quality, with intuition and stability in mind. The ATS Trinity™ imaging table exceeds the standards in patient positioning. The ATS Trinity™ surgical table features a C-arm compatible carbon fiber tabletop with a 550 LB lift capacity. The ATS Trinity™ is equipped with a motorized 4 way floating top to maximize patient imaging access. Other movements include Trendelenburg, Reverse Trendelenburg, Lateral Tilt, and Up/Down elevation.

Patient Safety

The ATS Trinity™ Surgical table is mobile on 4 anti-static castors. The floor lock system braces the table in place with unyielding stability. A battery backup system is a standard feature. The removable hand controller is ergonomic and outfitted with safety locks to avoid unintentional table movement. A return to level function is available as an option.

Movements

- 4 way floating top
- Left and right 15-degree lateral tilt
- Trendelenburg
- Reverse trendelenburg
- Up / down elevation


Features

- Carbon fiber floating top
- Silent electro-hydraulic movements
- 550 LB lift / weight capacity
- Electronic control safety lock
- Mobile anti-static casters with floor lock
- Floating top
- Battery backup
- Accessory rails
- Optional auto-reset / level function

INFINIUM®

Surgical Tables

Infinium ATS™ Surgical Tables



180°
Rotation


Slide

ATS Lateris™ Surgical Table

Features

- Silent electro-hydraulic movements
- 500 LB lift capacity
- Electronic control with safety lock
- Mobile anti-static casters with floor lock
- Horizontal slide for imaging applications
- 180-degree tabletop rotation
- Flex / reflex positioning
- Easy access to hydraulic system for service
- Optional battery backup
- Optional auto-reset / level function
- Optional 1000 lb lift capacity



180°
Rotation

ATS™ Surgical Table

Features

- Silent electro-hydraulic movements
- Electronic control with safety lock
- Mobile anti-static casters with floor lock
- Horizontal slide
- 180-degree rotation
- Manual / removable split leg section
- Easy access to hydraulic system for service

043.801.0002 REV1