

ClearVue®

Infinium ClearVue® VL3D

Disposable Blade Video Laryngoscope

The ClearVue® features a rechargeable lithium-ion battery for extended use. The HD screen offers a unique anti-fog and anti-shatter capability. The overall design is compact and ergonomic. A high-resolution view of the glottis will enable clinicians' first-attempt success while minimizing any chance of a complication during the intubation process.

COMPACT, ERGONOMIC, & DURABLE

- Lightweight - weighs less than 1 lb.
- Quick installation and removal of blade
- Quick shot camera button for recording
- 4 GB data storage (up to 80 hrs video)
- Optional HDMI output to remote display
- Screen is reinforced with anti-shatter protection
- Lightweight long lasting rechargeable battery
- 2-year warranty

COMPREHENSIVE BLADE OFFERING

- Macintosh blade sizes 1,2,3,4,5
- Macintosh channeled blade sizes 1,2,3,4,5
- More blades available

HIGH-DEFINITION

- 4-inch HD display
- Anti-fog lens
- 2.0 mega pixel camera with video



INFINIUM®

Video Laryngoscopes

Disposable Blade Video Laryngoscope

CLEARVUE™ TECHNICAL SPECIFICATIONS:

DISPLAY

Size	4"Full View
Monitor	LCD; Scratch-resistant and anti-glare screen
Resolution Ratio	640*480 (RGB (
Aspect Ratio	4:3
Video Refresh Rate	30FPS
Illumination	LED

LARYNGOSCOPE BLADE CAMERA INCLUDED

Resolution Ratio	2.0M Pixels
Field Angle	66°
illuminance	≥800LUX

BATTERY

Type	Rechargeable Lithium Battery
------	------------------------------

Long working time	≥120minutes
Voltage	3.7V
Capacity	3200mAh
Charging times	More than 300 times
Charging time	4 hours

POWER ADAPTER

Charging port	Micro USB
Input	100-250V,50/60Hz
Output	5V/ 2A

TRANSPORT/ STORAGE

Temperature	-10 (~ + 45 (
-------------	----------------

CONDITION

Humidity	≤93 (
AP	500hPa to 1060hPa

WORKING CONDITION

Temperature	5 (~ + 40 (
Humidity	30 (~85 (
AP	700hPa to 1060hpa

UNIQUE ANTI-FOG FUNCTION

Auto anti-fog function upon powering on

WEIGHT

321 g(disposable system); 376g (reusable system)
- light-weight, portable, comfortable

PHOTO/VIDEO FUNCTION

One click image capture 4G internal memory card
for video and images USB connection cable for
uploading images and video

HANDLE

ABS material - Comfortable ergonomic short
handle design

CLEARVUE ORDER INFORMATION & PART NUMBERS

00.300.0008	Laryngoscope 4in display kit with a set of disposable blades	306.300.0013	Laryngoscope Disposable Channeled Blade No.1
306.300.0008	Laryngoscope disposable blade No.1	306.300.0014	Laryngoscope Disposable Channeled Blade No.2
306.300.0009	Laryngoscope disposable blade No.2	306.300.0015	Laryngoscope Disposable Channeled Blade No.3
306.300.0010	Laryngoscope disposable blade No.3	306.300.0016	Laryngoscope Disposable Channeled Blade No.4
306.300.0011	Laryngoscope disposable blade No.4	306.300.0017	Laryngoscope Disposable Channeled Blade No.5
306.300.0012	Laryngoscope disposable blade No.5		

043.300.0003 REV2