

# Cleo™ Vitals & Capnography

## CLEO™ Vital Signs Monitor

- Blood pressure
- SpO2
- Touch screen
- Light weight
- Portable with battery backup
- Mount on rolling stand, IV pole, or wall

The Cleo™ is an intuitive approach to patient vital signs. It can be configured for any combination of non-invasive blood pressure and SpO2.

Weighing less than **3 LBS** the portable Cleo™ is well suited for any patient care area offering multiple vital sign combinations. It can act as a basic pulse oximeter/BP monitor or configured for NIBP/SpO2.

Cleo™ is well suited for both bedside and mobile spot check use. The monitor simplifies clinician use with a touch screen interface. A long-life lithium battery is standard & flexible mounting options are available.



Masimo SET® SpO2 technology available and recommended

### Exclusive NIBP features include:

**Inflation Pressure** – Users can set the initial inflation burst to a setting suitable for their patient base.

**Inflation Memory** – Recalls the patient's most recent Systolic measurement and adjusts the inflation pressure to the most comfortable level.

**Auto Inflation** – Users can set the Auto Inflation from 1-240mins, in 1min increments.

**Historical Data** – Recall last 300 readings.

**INFINIUM®**

**Patient Monitoring**

# Infinium Vital Signs & Capnography

## Cleo™ Capnograph

### Capnography Within Reach

- EtCO<sub>2</sub>
- Optional SpO<sub>2</sub>
- Optional blood pressure
- Touch screen
- Light weight
- Portable with battery backup
- Mount on rolling stand, IV Pole, or wall

### Flexible Configurations

- EtCO<sub>2</sub> - mainstream or sidestream
- Blood pressure
- SpO<sub>2</sub>

Weighing in at less than 3 LBS the portable Cleo™ is well suited for any patient care area by offering flexible vital sign combinations (EtCO<sub>2</sub>, BP, SPO<sub>2</sub>).

The Cleo™ simplifies clinician use by incorporating a touch screen with a simple user interface making the Cleo® intuitive for any user. A long-life lithium-Ion battery is standard and many mounting solutions are available for the Cleo™.



*"Reliable, accurate, and the lowest cost per patient capnograph."*



**Pain Management**



**Procedure Sedation**



**Sedation Dentistry**

**INFINIUM®**

043.500.0003 REV1

Infinium Medical  
12151 62nd St North #5  
Largo, FL 33773 • USA

☎ (1) 866-918-8434  
🌐 [infiniummedical.com](http://infiniummedical.com)  
✉ [sale@infiniummedical.com](mailto:sale@infiniummedical.com)