

INFINIUM[®] FOCUS[™]

Powerful Therapy. Precise Results.



*Deep Penetrating Focused
Hybrid Electromagnetic Shockwaves*

Therapy Equipment

Infinium FOCUS™

Hybrid Shockwave Therapy System

The Infinium FOCUS™ Focused Adjustable Shockwave Therapy System is designed to deliver deep penetrating focused shockwaves with unmatched precision. Engineered for versatility and results, it supports a wide range of therapy types—from musculoskeletal rehabilitation to advanced pain management.



KEY FEATURES

- » Clinical Grade Strength
- » Energy Flux Density 0.06 – 0.39 MJ/MM2
- » 29 Pre-Programmed Therapy Protocols
- » Easy to Operate Touch Screen User Interface
- » Adjustable Depth Focused Shockwave (5 - 60mm)
- » Pain Free Quick Treatments

INTUITIVE TOUCH SCREEN



INFINIUM®

Infinium FOCUS™ - ADVANTAGES



ADJUSTABLE DEPTH PENETRATION

Focused shockwaves can reach tissue depths from .5 to 6cm with overall penetration up to 12cm.



PRECISE TARGETING

Energy can be concentrated at a specific focal point, ideal for treating localized injuries or chronic pain.



HIGH ENERGY LEVELS

Capable of delivering stronger therapeutic energy (MJ/MM²) for stubborn or deep-seated conditions.



VERSATILITY IN THERAPY

Effective for musculoskeletal, orthopedic, neurological, urological applications and more.



FASTER RESULTS

Patients often experience quicker pain relief and accelerated healing compared to radial treatments.



REDUCED DISCOMFORT

Capable of delivering stronger therapeutic energy (MJ/MM²) for stubborn or deep-seated conditions.



SCIENTIFIC BACKING

Focused shockwave therapy has a longer history of clinical research and proven outcomes across multiple medical specialties.



ELECTRO-HYBRID Shockwave Therapy

FOCUS™ Shockwave Therapy System Specifications

System Type	Hybrid Electromagnetic focused pulse therapy
Treatment Depth	5m to 120mm
Focal Depth	5mm to 65mm
Frequency	1-10 Hz
Energy Flux Density	0.06 – 0.39 MJ/MM2
Packaging	Aluminum Carrying Case
Protocols	29 preset treatment recommendations
Controls	Color touch screen for all software operations
Applicators	Water filled bulb with controls
Applicator Lifecycle	500,000 shocks per handpiece 1,000,000 shocks per pulse
Dimensions	23.5 × 18 × 8 in (L/W/H)
Unit Weight	31 kg / 68 lbs
Warranty	2 yrs controller unit and 500,000 pulses for applicator

25+ PRESET MEDICAL PROTOCOLS

- Chronic Back Pain
- Calcific Tendinitis of the Shoulder
- Degenerative Joint Ankle
- Patellar Tendinopathy
- Tendonitis of the Wrist
- Tennis Elbow
- Hip Pain
- Tibial Stress Syndrome
- Heel Spurs
- Erectile Dysfunction
- Trigger Points Neck
- Upper Trapezius Strain / Spasm
- Cervical Spondylosis
- Bicep Tendonitis
- Rotator Cuff
- Golfer's Elbow
- Carpal Tunnel Syndrome
- De Quervain's Tenosynovitis
- Lumbar Spondylosis
- Herniated Nucleus Pulposus
- Paravertebral Syndrome
- Lumbar Mobilization Block
- Piriformis Syndrome
- Hamstring Shortening
- Osteoarthritis of the Knee
- Calf Muscle Strain
- Quadriceps Spasm



Infinium Medical
12151 62nd St North #5
Largo, FL 33773 • USA

- ☎ (1) 866-918-8434
- 🌐 infiniummedical.com
- ✉ sale@infiniummedical.com

