

HS360 SurgiStool

INFINIUM[®]



Ergonomic Hydraulic Surgical Stool

HS360 SurgiStool

Unmatched in Ergonomics Function & Efficiency

2 YEAR
WARRANTY

Infinium's HS360 surgical stool is an exceptionally versatile and reliable seating solution, perfectly suited for a wide range of surgical procedures. Its innovative design incorporates numerous features and benefits to enhance user comfort and efficiency. Equipped with intuitive hydraulic foot pedals for effortless height adjustments, along with easily adjustable backrest and armrests, this SurgiStool ensures optimal ergonomics and seamless operation during demanding tasks.

- Adjustable contoured backrest for personalized support
- Dual-density foam seat cushion for enhanced comfort
- Foot-controlled hydraulic system with Y-type base and four smooth-rolling casters for easy mobility
- Front-end locking system for secure positioning
- 360° rotating armrest pads with 180° armrest support for flexible adjustments



| | |
|---------------|-------------------------------------|
| High Position | 28 in \pm 1 in (71 cm \pm 2 cm) |
|---------------|-------------------------------------|

| | |
|--------------|-------------------------------------|
| Low Position | 21 in \pm 1 in (53 cm \pm 2 cm) |
|--------------|-------------------------------------|

| | |
|-----------------|------------------|
| Weight Capacity | 300 lbs (136 kg) |
|-----------------|------------------|

Smooth Hydraulics



Three Seat Shapes



INFINIUM[®]

Infinium Medical
12151 62nd St North #5
Largo, FL 33773 • USA

 (1) 866-918-8434
 infiniummedical.com
 sale@infiniummedical.com