

**LeMans 300 SERIES**

**INFINIUM<sup>®</sup>**



**Barrier-Free Examination Tables**

# LeMans 300

## Unmatched in Ergonomics Function & Efficiency

Developed from years of research and feedback from hundreds of healthcare professionals, we have designed the LeMans Series to be the most usable and patient-friendly table in its class.

The LeMans 300 is **compact** but still offers features normally found only on more expensive tables with larger footprints.

- » Load capacity of 450 pounds
- » Ideal height adjustment from 19" to 39"
- » Integrated conforming leg rest
- » Seamless easy to clean upholstery options
- » Integrated retractable stirrups
- » Pass through rear drawer
- » Front storage drawer
- » Position memory settings



20 Special Order Upholstery Colors



Foot & hand control of  
height and back rest.




19" to 39"



Up to 84" Flat




# Infinium LeMans 300 Options & Accessories



ADJUSTABLE FLUSH  
LEG REST




PAIR OF RAILS



WARMING DRAWER



WIRELESS PEDAL



STAINLESS STEEL TRAY  
(REPLACING THE STD TRAY)



PAIR OF STRAPS  
FOR GOEPPLE LEG  
SUPPORTS



INTEGRATED MAINS SOCKET



SLIDING RAIL VICE



ROTARY RAIL VICE




ARM SUPPORT




ARM SUPPORT COVER



PAIR OF STIRRUPS



ADJUSTABLE ROLL  
HOLDER



HEADREST



SELF-TIGHTENING DRIP  
STAND W/ 2 HOOKS



DRIP STAND  
W/ 2 HOOKS



LAMP SUPPORT LID  
FOR RAIL



LAMP HOLDER



RVS100  
LEFT BRACKET



RVS100  
RIGHT BRACKET



# LeMans 300 Exam & Procedure Table Technical Specifications

**Load Capacity** 450 lbs

**Height adjustment** 19" to 39"

**Net Weight**

Unit: 412 lbs

Upholstery: 26.5 lbs

**Gross Weight**


Unit: 485 lbs

Upholstery: 35 lbs

**Package Volume**

Unit: 24.7 cubic ft

Upholstery: 5.7 cubic ft



## INFINIUM®

Infinium Medical  
12151 62nd St North #5  
Largo, FL 33773 • USA

(1) 866-918-8434  
infiniummedical.com  
sale@infiniummedical.com